

Dual Lumen Catheter



This can be used to inject contrast, insert guide wire and can also be used as a dilator. Can be used for a number of Functions such as maintaing drainage, infusion of contrast media and ureteral dilation.



Key Features :

- Provides an independent guide wire lumen and injection lumen
- Tapered tip enables the dilation of ureter
- Smooth surface for ease of introduction
- Eliminates the need for multiple catheterization

Ordering information:

Size	8 fr.	10 fr.	12 fr.
Length			
50 cm	BUDC0108050	BUDC0110050	BUDC0112050