

INDOSCOPE

bioradmedisys™
science for people

e-lite

Features

- Compact video system with 15.5 inch integrated screen
- Integrated optical connector
- Built in Light Source
- Multiple Video Output ports
- Records Videos and Photo on USB drive



Technical Specification

- Storage / Data Transfer - USB 3.0
- Image Functions - White Balance, Brightness Control, Image Capture, Video Capture
- Video Connections - DVI , SDI & HDMI
- Power - 100 V AC – 240 V AC
- Lamp - 30 W LED
- Image Resolution - 1920 x 1080 HD

Indoscope Processors



Features	Control Unit	CVS e-lite	CVS
Built In Light Source	●	●	●
Video Recording	●	●	●
Screen	—	15.5 Inch	21 Inch
BNB	●	—	—
BCR	●	—	●
Optical Connector	External	Integrated	External
BNB- Blue Narrow Band BCR- Blood Colour Reduction			
Connects To			
Indoscope Sleek (7.5 Fr)	●	●	●
Indoscope ACE (5.5/7 Fr)	●	●	●
Indocam	●	—	●
Video Output			
DVI	●	●	●
HD-SDI	●	●	●
HDMI	●	●	—

Warranty : Biorad Medisys shall warrant the products against defective material and faulty workmanship for a period of one (1) year for control unit, cvs, e-lite except lamps of Endoscopic light sources, rubber parts, protection sheaths, glass parts, cables, screens, light guide cables from the date of shipment of these products

The warranty shall not be applied to a) improper use, maintenance and/or negligence by customer, and b) improper repair, adjustment or negligence by any other party not authorized by BIORAD MEDISYS.

Disclaimer :

Note: Specifications indicated in this brochure are subject to change without any prior notice.

bioradmedisys™
science for people

Corporate Office & Manufacturing Unit
Plot No.48/3 & 48/7, Pashan-Sus Road, Sus Village,
Near Tapkir Vasti, Pune - 411 021. Maharashtra, India.
Phone: 020-35000200

www.bioradmedisys.com

E: [biorad@bioradmedisys.com](mailto: biorad@bioradmedisys.com)

Follow us on [f](#) [t](#) [in](#) [v](#)