





TOMI
MORCELLATOR



EFFICIENT & CONVENIENT PROSTATE MORCELLATION





Technical Specifications

Blade Diameter	4.5 mm
Blade Shaft Length	400mm
Max. Morcellation Speed	1200 rpm
Max. Output Torque	6.0 mN m
Suction	Integrated
Vacuum Pressure	0.08 Mpa
Water Flow Rate	1500 ml/min
HandPiece Sterilisation	Autoclave
Foot Switch	Two Stage
Accessories	Aspiration tube set, Waste container, Blade
Power Supply	230V~ 50Hz/60Hz

Warranty : Biorad Medisys shall warrant the products against defective material and faulty workmanship for a period of one (1) year except rubber parts, protection sheaths, glass parts, cables, from the date of shipment of these products

The warranty shall not be applied to

- a) improper use, maintenance and/or negligence by customer, and
- b) improper repair, adjustment or negligence by any other party not authorized by BIORAD MEDISYS.

Disclaimer :

Note: Specifications indicated in this brochure are subject to change without any prior notice.



Corporate Office & Manufacturing Unit
Plot No.48/3 & 48/7, Pashan-Sus Road, Sus Village,
Near Tapkir Vasti, Pune - 411 021. Maharashtra, India.
Phone: 020-35000200

www.bioradmedisys.com

E: biorad@bioradmedisys.com

Follow us on    